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DISPOSAL PROCEDURE (Dir. 2012/19/Ue-Weee) The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. **Caution!** The wrong disposal of electric and electronic equipment may involve sanctions.

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
GB

1. 2000 ml Vessel
2. Vessel plug
3. Protection device
4. ON-OFF button
5. Vacuum regulator
6. Vacuum gauge
7. Multi-voltage power supply unit and battery charger socket
8. 20 cm silicon tube
- 8a. 10 cm silicon tube
9. Antibacterial filter
10. Working light
11. Light ON-OFF button
12. Accessory compartment
13. Ambulance mounting bracket

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

GB

1. **Flashing red Led:** the battery is no longer capable of supplying maximum charge duration; replace it.
2. **Red Led permanently on:** the battery is flat, use the appliance for a maximum of 3 minutes and re-charge the battery
3. **Yellow Led permanently on:** the battery is in reserve (about 10 minutes' operating time left); recharge it if possible.
4. **Green Led permanently on:** appliance operating, battery charged.
Green Led flashing: battery under charge.

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
Fig. - Abb. 7e


Fig. - Abb. 7f

The MINIASPEED BATTERY EVO PLUS aspirator is a professional portable device suitable for use in ambulances, specific for secretion removal. It has a vacuum regulator, a vacuum gauge and 2000 ml container with a protection device to prevent liquid from entering the suction pump, which cuts off the aspiration flow. It has lifetime lubrication and is easy to handle, simple to use, reliable, strong and silent. This model has a 12 V rechargeable battery; when connected to the ambulance mounting bracket it not only operates on the ambulance's battery but also recharges. The MINIASPEED BATTERY EVO PLUS aspirator has the following 3A accessories: 2000 ml container with protection device, power supply cord with multi-voltage power supply unit, 10 cm long sterilisable silicone connection tube, 20 cm long sterilisable silicone connection tube, 130 cm long sterilisable silicone connection tube, **sterile, disposable** cannula, **sterile, disposable** manual flow regulator and **disposable** antibacterial filter.

N.B.: Only use original 3A accessories.

IMPORTANT WARNINGS

This is a medical device and must be used by qualified staff. It must be operated as indicated in this user instruction manual. It is important for the operator to read and understand the information on use and maintenance of the unit. If you have any questions, contact your stockist. MICROBIAL CONTAMINATION: in the case of illnesses with a risk of infection or microbial contamination, the accessories should be thoroughly cleaned and sterilised after each use.

The manufacturer has made every effort to ensure that all its products are of the highest quality and are reliable. Nevertheless, as for all electrical appliances, fundamental safety standards must be observed:

- Children and persons who are not self-sufficient may only use the unit under the strict supervision of a responsible adult who has read this manual.
- The device must always be used by specifically trained staff who have read this manual.
- Never use adapters for supply voltages different to the voltage shown on the data plate on the back of the unit. Keep the cord away from hot surfaces.
-  The unit should not be used in the presence of inflammable anaesthetic mixtures with air, oxygen or nitrous oxide.
- Never handle the supply cord plug with wet hands or use the unit when taking a bath or a shower. Never leave the unit close to water, do not submerge it in any liquids. If by chance it falls into water, pull the plug out of the socket immediately before recovering it. Do not use the unit if the plug or supply cord are worn or wet (send it immediately to your stockist).
- Although the device is protected against sprays of water, do not use it in the rain.
- Only authorised personnel may perform maintenance and/or repair work. Unauthorised repairs annul the guarantee.
- Ensure that the connections and container closing are performed carefully to prevent suction losses.
- Do not tip the container over while connected to the unit when operating, as liquid may be sucked inside the appliance causing pump damage. If this happens, switch the aspirator off immediately and empty and clean the container. Send the appliance to your stockist.
- When the protection device intervenes aspiration is interrupted; empty the container and perform the cleaning operations.
- The cannula and the manual control of the aspirated flow are **sterile, disposable** products: They must be replaced after every application.
- Check the use-by-date on the original packaging of the cannula and check the integrity of the sterile packaging. If expired and/or deteriorated, replace it.
- The **disposable** antibacterial filter must be replaced after every application.
- Never use the battery charger with other appliances or for uses other than that established by this manual. Never use MINIASPEED BATTERY EVO PLUS with other power supply units.
- In view of their length, the power cord and connection hose could constitute a strangulation hazard.
- Use of the device in ambient conditions other than those specified in the manual may seriously impair its safety and technical characteristics.
- In the event of aspiration without the container and/or antibacterial filter, or if it is suspected that substances have entered the aspiration circuit, contact your stockist at once.
- Never modify the device's mounting bracket or the bracket installed in the ambulance as this may impair the safety of both the device and its installation environment.
- Always use the device in the vertical position on an unobstructed, stable, flat surface.
- The lifetime of the device itself is 5 years, and that of the container and silicone tubes 1 year or 30 sterilisation cycles. However, the lifetime may vary depending on the operating environment.

USE INSTRUCTIONS

Each time before use, ensure that all the accessories are perfectly clean according to the instructions indicated in the "CLEANING AND DISINFECTION OPERATIONS".

1. OPERATION WITH MULTI-VOLTAGE POWER SUPPLY UNIT

- 1.1 Check the charge level of the ambulance's battery before using the device
- 1.2 For use outside the ambulance (i.e. not fixed to the bracket [19]) with the internal battery flat or not full charged, connect the multi-voltage power supply (16) to the socket on the device (7) and connect the cord (17) to a mains power socket.
- 1.3 Connect the appliance as indicated in Figure 1 (page 1).
- 1.4 Start the appliance by placing the switch to the "I" (ON) position (4) (green pilot light on).
- 1.5 The vacuum regulator (5) can be used to set the level of vacuum required (bar/KPa). Turn the knob clockwise, in the "+" direction, to obtain a higher vacuum or anticlockwise, in the "-" direction, for a lower vacuum; these values are read on the vacuum gauge (6).
Important: the vacuum values on the control decal are purely for guidance; always refer to the vacuum gauge reading.
- 1.6 After the application, switch off the device, remove the supply cord (17) from the mains socket and disconnect the power supply unit (16) from the device. Perform the cleaning operations as described in the "CLEANING AND DISINFECTION" paragraph.

2. OPERATION WITH BATTERY AND BATTERY CHARGER WITH MULTI-VOLTAGE POWER SUPPLY UNIT

- 2.1 The device is supplied with the battery partially charged. We therefore recommend charging it before use.
- 2.2 To charge the battery, with the device switched off connect the multi-voltage power supply unit (16) to the device's socket (7) and to the mains power supply by means of the cord (17). Charging time: about 4 hours. Charge duration: about 40/45 minutes at maximum suction power.
- 2.3 Operation with internal battery only:
Switch the device on by pressing the "ON-OFF button (4) (green light comes on). If the green light goes out during use, and the yellow light comes on, there is about 10/15 minutes of battery charge left (reserve level). Therefore: terminate the application if possible. If the application cannot be terminated, it can be continued until the red (battery flat) light comes on, but if this occurs do not use the appliance for more than 3 minutes to avoid damaging the battery. To continue the application with the battery flat (red light on), connect the multi-voltage power supply unit (16) (as described in point 2.2).
- 2.4 At the end of the application, switch the device off by pressing the button (4). Connect the power supply unit (16) to charge the battery (as described in point 2.2). If the device is not in use and is not connected to the bracket (19) in the ambulance, the power supply unit should be left connected (16) to ensure that the battery is always charged to the optimal level.
- 2.5 For aspiration procedures see points 1.3; 1.4; 1.5.

N.B.: when the battery power drops below a set level the aspirator shuts down to avoid damage to the battery. If absolutely necessary, the user can restart the aspirator for 1 minute by pressing the ON/OFF button again.

3. OPERATION WITH THE DEVICE INSTALLED IN THE AMBULANCE

- 3.1 The device can also operate connected to the ambulance mounting bracket (19).
- 3.2 With the device connected to the bracket (19), there is no limit to patient treatment time because the device is powered directly by the ambulance battery.
- 3.3 The device's battery is charged by leaving it connected to the ambulance mounting bracket (19).
 - 3.3.1 The device's battery is also charged while treatment is in progress (with the device connected to the ambulance mounting bracket).
- 3.4. For aspiration procedures see points 1.3; 1.4; 1.5.

4. SECRETION VESSEL - 2000 ml (Fig. 3)

The secretion vessel set consists of an overflow valve, a vessel (1) in clear material (polycarbonate) and a blue plug (2).

Connection: connect one end of the 10 cm tube (8a) to "VACUUM" connector of the blue top (2) and the other end connect it to the antibacterial filter connector (7). Then connect one end of the 20 cm tube (8) to the free connector of the antibacterial filter (7) and the other end to the "INLET" connector of the aspirator. The antibacterial filter also protects the aspiration circuit from any contaminating agents sucked in during use. Do not use the aspirator without the antibacterial filter, because from a bacteriological point of view, it becomes dangerous for the patient. Keep the device vertical to allow the overflow to function correctly. Connect one end of the 130 cm tube (9) to the "PATIENT" connector of the blue cap (2), while at the other end to connect the disposable, sterile manual regulator (11) and connect the disposable, sterile cannula (10) to this latter.

All the components of the vessel can be sterilised using a conventional system in an autoclave at a temperature of 121°C, or by boiling for 10 minutes. We recommend replacing the complete vessel at every 30 sterilisation cycles. Do not overturn the vessel during use, in order to prevent the intervention of the non-return valve (3-4-5); should this occur, switch the aspirator off and detach the tube connected to the antibacterial filter. Never use the aspirator without the secretion collection vessel and/or without the antibacterial filter.

DESCRIPTION OF USE OF WORKING LIGHT

The MINIASPEED BATTERY EVO PLUS is fitted with a LED working light (10) which allows use of the device in dark conditions and/or at night.

Press the button (11) to switch on the light, as shown in Figure 5 of page 3; press the button (11) again to switch it off. The light switches off automatically after about 60 minutes.

The light is of LED type, so the amount of current it consumes is minimal and it can be used even during treatment; however, it should be switched off whenever it is not required and/or at the end of the treatment.

CLEANING AND DISINFECTION OPERATIONS

N.B.: If using chemical disinfectants, follow the manufacturer's instructions exactly.

- The cannula and the aspirated flow manual command are **sterile, disposable** products and must be replaced after every application.
- The antibacterial disposable filter must be replaced after every application.
- Never leave the appliance in water or submerged; clean the external casing of the appliance using only a damp cloth with detergent (non abrasive).

DESCRIPTION OF AMBULANCE MOUNTING / RELEASE SYSTEM

1. Fit the bottom of the bracket fixed to the aspirator into the recess provided in the bottom of the bracket mounted in the ambulance (as shown in Figure 6 of page 3).
2. Keeping the device engaged with the bracket, push it until it engages with the hook in the top of the yellow mobile bracket (as shown in Figure 6a):

N.B.: to ensure that the device engages correctly, press the mobile bracket, push the device fully down and release the mobile bracket (as shown in Figure 6b).

 **ALWAYS CHECK THAT THE DEVICE HAS ENGAGED WITH THE BRACKET.**

3. To release the device, press the yellow mobile bracket and lift the aspirator up and towards you (as shown in Figure 6c).

 **KEEP A FIRM GRIP ON THE HANDLE OF THE ASPIRATOR WHILE RELEASING IT.**

DESCRIPTION OF PROCEDURE FOR MOUNTING BRACKET IN AMBULANCE

The ambulance mounting bracket is supplied assembled and complete (as shown in Figure 7):

1. Wall mounting bracket.
2. Mobile bracket.
3. Mobile bracket spring.
4. Block complete with electrical contacts and cord.
5. 2 metre power cord for connecting the ambulance battery.

Mounting the bracket in the ambulance

- The bracket must be installed on the panel inside the ambulance by skilled staff authorised by the vehicle's constructor/outfitter.
- The fixing points are shown in Figure 7a, Figure 7b and Figure 7c.
- M5 Allen screws (UNI 5931-DIN 912-ISO 4762) should be used since they undergo deceleration testing at 10g as required by the EN 1789 standard. Otherwise, screws (with the same mechanical characteristics) having (non-tapered) heads, max. Ø10, and no more than 7 mm in height (including the washer, if fitted) can be used.

- See Figure 7d for the fixing centre distances.

Wiring connection/polarity

- The bracket is supplied complete with the ambulance battery connection cord (5) (already fitted), which powers the device and charges its battery when it is connected to the bracket.
- Connect the wires to the ambulance following the polarity shown in Figure 7e.
- If the cord is disconnected from the bracket for any reason, reconnect it following the polarity shown in Figure 7f.

N.B.: Take great care over the polarity (+ and –) because if it is reversed the device's electronic circuit board and battery will be damaged beyond repair.

PROBLEMS, CAUSES AND SOLUTIONS

| PROBLEMS | POSSIBLE CAUSES | SOLUTIONS |
|--|---|--|
| Excessive noise. | Damaged pump or blockages in the internal aspiration circuit. | Send to the assistance circuit. |
| The unit switches on but does not aspirate. | - Damaged pump. - Vacuum regulator fully open. Connection tubes disconnected and/or badly connected, broken connection tubes. Container not in a vertical position, full, or defective overflow valve. Possible blockage of the hydraulic circuit inside the unit. | - Send to the assistance circuit. - Check the position of the vacuum regulator. Check the connections and the integrity of the tubes. Position the container in a vertical position, check the overflow valve (blocked) and/or replace the silicon tubes. |
| The vacuum rate cannot be regulated. | Damage to the internal hydraulic system or blockage of the connection tubes to the aspiration unit. | Send to the assistance circuit. |
| When the appliance is switched on, the protection fuse always trips. | Pump probably damaged or in shortcircuit. | Send to the assistance circuit. |
| The vacuum gauge does not work. | Liquid penetrating the pneumatic circuit. | Send to the assistance circuit. |

Note: if you experience faults or malfunctioning problems different to those listed above, always and exclusively contact authorised assistance centres.

TECHNICAL CHARACTERISTICS

Professional portable medical device for secretion removal, suitable for fixed and transportable use. Complete with flame retardant ABS V0 casing and piston-cylinder electrical compressor lubricated for life.

Risk class under Directive 93/42/EEC:

Suction class:

Protection level against entry of liquids and solids:

Voltage:

Absorption:

Adjustable vacuum level:

Vacuum gauge precision class:

Air flow:

Temporary use:

Dimensions:

Weight:

Noise level:

Multi-voltage switching battery charger:

Internal battery:

Battery charge duration:

Operating conditions:

Temperature:

Air humidity:

Storage conditions:

Temperature:

Air humidity:

Operating-storage conditions:

IIb

High vacuum / High flow

IP44D

12VDC

3A

from 0 to –0.85 bar

2.5% under UNI EN 837

30 lt/min operating without restrictions +/- 10 %

max. 60 minutes

38 (length) x 13.5 (width) x 35 (Height) cm

4.9 kg approx

61dBA

PRI: 100 ÷ 240V~ 50 – 60Hz

SEC: 14VDC 4.28A

12VDC 4Ah Lead, hermetic

45 minutes at maximum suction power

min. 0° C; max 40° C

min. 10 %; max 95 %

min. –10° C; max 50° C

min. 10 %; max 95 %

min. 690 hPa; max 1060 hPa

SYMBOLS USED

Type BF appliance

Attention, check the use instructions

Switch On

Switch off (or battery on charge)

Direct current

Alternate current

Do not use the unit when taking a bath or a shower

CE 0434 Conforms to Directive 93/42/EEC for medical devices

The device contains a hermetic lead battery. It must be disposed of in accordance with current regulations on the disposal of toxic-harmful waste

IP44D

The device's casing is protected against solid particles having diameter of 1 mm or above, spraying water and access to hazardous parts with a wire

Thermally protected safety transformer

Sterilisation by ethyl oxide

Disposable

Electromagnetic compatibility

additional notes in conformity with the IEC 60601-1-2/A1:2004-09 Standard

| EMISSION | | |
|--|--------------------|--|
| Emission test | Conformity | Electromagnetic Environment - guidance |
| Emissioni RF Cispur 11 | Group 1 | MINIASPEED BATTERY EVO PLUS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| Emissioni RF Cispur 11 | Class B | MINIASPEED BATTERY EVO PLUS is suitable for use in all establishments and those directly connected to the public low-voltage power supply network that supplies buildings used to domestic purposes. It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, adopting provisions in the installation, such as a longer distance from potentially sensitive appliances. |
| Harmonic Emissions IEC 61000-3-2 | Class A Conforming | It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings user for domestic purposes. |
| Voltage fluctuations/ flicker emissions IEC61000-3-3 | Conforming | |

| IMMUNITY ASPECTS | | | |
|--|--|--|---|
| MINIASPEED BATTERY EVO PLUS is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | | |
| Immunity test | Test level EN 60601-1-2 | Compliance Level | Electromagnetic environment - guide |
| Electrostatic discharge (ESD) EN 61000-4-2 | ± 6kV contact ± 8kV air | ± 6kV contact ± 8kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Burst/Fast Transient EN 61000-4-4 | ± 2kV power supply lines | ± 2kV power supply lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge EN 61000-4-5 | ± 1kV differential mode ± 2kV common mode | ± 1kV differential mode ± 2kV common mode | Mains power quality should be that of a typical commercial or hospital environment. |

| | | | |
|--|--|--|--|
| Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11 | < 5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds | < 5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds | Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions. It is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency magnetic field EN 61000-4-8 | 3 A/m | 3 A/m | Magnetic power frequency fields should be that of a typical commercial or hospital environment. |

IMMUNITY ASPECTS AT R.F.

MINIASPEED BATTERY EVO PLUS is intended for use in the electromagnetic environment specified below. The customer or the user of the navigator should assure that it is used in such an electromagnetic environment.

| Immunity test | Test level EN 60601-1-2 | Compliance Level | Electromagnetic environment - guide |
|---------------------------|-----------------------------|-----------------------------|--|
| RF conducted EN 61000-4-6 | 3 Veff from 150kHz to 80MHz | 3 Veff from 150kHz to 80MHz | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from that equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ from 150kHz to 80MHz $d = 1.2 \sqrt{P}$ from 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ from 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). |
| RF conducted EN 61000-4-3 | 3 V/m from 80MHz to 2.5GHz | 3 V/m from 80MHz to 2.5GHz | |

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE DEVICE

MINIASPEED BATTERY EVO PLUS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of the transmitter (m) | | |
|---|---|--|---|
| | From 150kHz to 80MHz $d = 1.2 \sqrt{P}$ | From 80MHz to 800MHz $d = 1.2 \sqrt{P}$ | From 800MHz to 2GHz $d = 1.2 \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Notes:

- (1) At 80 MHz and 800 MHz the separation distance for the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.