

Massage system Pro6-MK400L

For professionals



Features:

- The ideal massage system for treatment of clients with lymphedema, injuries, etc.
- Comes with two leg cuffs, size XL.
- Uses compression therapy with 6 chambers.
- One channel model for two cuffs.
- Extensive adjustment possibilities.
- Each chamber can be switched on/off separately.
- Improves lymph activity, blood circulation and overall condition.
- Various additional cuffs available.
- Not labor intensive.
- Interesting price, fast payback times.



How it works:

The system stimulates the lymphatic system and blood circulation system by massaging legs, arm, hips and abdomen with air-filled cuffs. Waste products in the muscles are removed quickly. The massage can treat fat deposits, hardening of veins, skin and muscle tissue. It helps against cellulite and gives a nicer skin. And above all: it provides a deep relaxation.

What is DoctorLife:

DoctorLife is a registered trademark of Daesung Maref Corp. Ltd. from South Korea. They make air-driven massage (IPC) equipment for medical treatments, sports massage and wellness applications. The systems are sold in more than 70 countries. North Sea Health Care BV is an official importer.

More info:

For more information please contact:

North Sea Health Care BV

Maarnse Grindweg 22

3951LK Maarn

The Netherlands

Tel: +31 (0)6 22906310

E-mail: info@doctorlife.nl

Visit our website: www.doctorlife.nl