

FlexiHealth, personal protective screen for medical use

Reinforce safety with your patient



Description: Personal protective equipment to reduce the risk of viral exposure in medical procedures. Intended for dentists, **ADDITIONAL ADVANTAGES** ophthalmologists, plastic surgeons, etc.

FlexiHealth acts as a barrier to protect against droplets and aerosols that cause transmission of infection. It reinforces the safety of the health professional without altering the working conditions with the patient.



Characteristics: It consists of a 3-4 mm thick polycarbonate screen, which provides optimal visibility of the operating field, as well as extreme resistance to shocks and scratches.

It is equipped with a flexible articulated arm and an adjustable clamp that allows it to be fixed to any surface, according to the user's needs, occupying a minimum of space.

ADDITIONAL ADVANTAGES



Specifications:

Arm length: 1000 mm

Screen sizes:

410 x 410 mm

330 x 330 mm