

Work happily



V1.0


naggura'swop

naggura.com

Our products



S703LB S703DB




The Swop 703 is an electrical 3-section couch with a one-piece mattress and wooden structure.

Wood, Metal or Pmma




Collections


Nordic
METAL



Yellow




Blue




Pink

Foll
METAL




Ash




Clay

Iron
METAL



Inox



Black

Om
WOOD



Oak




Wenge

Techno
PMMA



Fluor




Fuxia



Ice

S703






S703LB



S703DB

Dimensions



General Specifications

MOTORS 3 engine (height, footrest, backrest)
DIMENSIONS (LENGTH X WIDTH) 200 x 75 cm
HEIGHT (MIN. AND MAX.) 59/85 cm
WEIGHT 110 kg

THICKNESS 10 cm
DENSITY 30 kg/m ³
LIFTING CAPACITY 175 kg
OPERATION Hand controller (included)/ Foot controller (optional)

POWER SUPPLY 230 V - 50/60 Hz
MAX. INTENSITY 2 A
ELECTRIC MOTOR 6.000 N

CUSHION Upholstery Seamless round edges 100% Vinyl-compound Permaguard 3 advanced vinyl protection (against germs and stains), UV re- sistant, superior abrasion 100K
--