

EQUIPAMIENTO OPCIONAL

1. Drip bag support
ref. MB-DIAMED-01.
2. Table ref. MB-DIAMED-02.
2. Adjustable head rest.
ref. MB-DIAMED-03.
2. Push bar
ref. MB-DIAMED-04.



ELECTRICAL CHARACTERISTICS

Electric feed	230V 50Hz
Maximum consumption	1.6 A, 370 Volt
Protection Indicator	IP 66
Class of protection	Class I
Level of protection in event of surge	Type B

TECHNICAL DATA

External dimension	183.7 x 90.6 cm
Back rest dimensions	59.6 x 817 cm
Pelvic area dimensions	59.6 x 55.5 cm
Leg rest dimensions	47.2 x 59.6 cm
Back rest inclination	17°
Trendelenburg dimensions	13°
Extendible foot rest	42.8 - 57.8 cm
Maximum load	135 kg
Weight excluding accessories	45 kg



NOTE: The dimensions and characteristics of the products presented may be changed without notice.



Hospital Furniture DIAMED multifaceted chair

www.medical-iberica.com

Lanzahita, 6 - Pol. Ind. Albarreja - 28946 Fuenlabrada, Madrid (Spain) - T. 91 208 94 00 - F. 91 208 94 01

DIAMED CHAIR

The multifaceted **DIAMED** chair has been developed to meet the great demands of hospitals around the world.

The attractive design of the **DIAMED** chair, helps create a feeling of wellbeing for the patient. The exclusive design establishes the highest of standards, thus increasing the prestige of any hospital or medical centre.

BRIEF OVERVIEW

The design allows for easy decontamination and disinfection.

The pleasing design is combined with compact proportions.

Makes for a comfortable place to relax and thus, a more pleasurable experience for the patient.

Reduces the physical effort required to move the patient.

EASY CLEANING

Our products are manufactured with materials of the highest quality, non porous to guarantee long term resistance and minimise risk of infections and reduced the cleaning time.



EASY MOVEMENT CONTROLS

The chair's control panel has clear pictograms for the patient and medical staff. Not only does the control adjust the chair height, the incline of the back rest and the leg section, it also includes two pre-programmed functions which allow the chair to be moved into cardiac seat or Trendelenburg positions at the touch of a button.



IMPROVING THE PATIENT'S POSITION AND RECOVERY

The wide range of positions that the seat has, allows the patient to find the most comfortable position for their own needs.



Foot rest extension.

EASE OF MOBILITY

Its four castors of 125 mm in diameter facilitate movement, reducing the physical effort required to move it.



Trendelenburg 16°.

DIAMED

A chair that adapts to each user's individual needs



ARM RESTS

The adjustable ergonomic arm rests facilitate access. The DIAMED chair system has been developed to allow both the extension of the arm rests and multiple positioning.

COMFORT

The **DIAMED** chair has undergone extensive testing and research in laboratories to achieve the maximum ergonomic adaptability and optimum comfort, necessary to improve the quality of life and health of patients.

100% of surface area is usable.

Base planks are jointed by way of low tension electric motors in accordance with the guideline EN 60601-1.

UPHOLSTERY

Seat and backrest in high density PU foam. Upholstered in autoextinguishable vinyl fabric, in accordance with regulation M1. Easy to clean and maintain.



HANDLING

To facilitate movement, the chair has four 125 mm castors, a centralised brake system and power steering.

BASE

Frame constructed using steel tubing and coated in epoxy paint.

Elevation system which provides great stability and a load of up to 230 kg.

Incorporates supports which allow for additional accessories.