

## TECHNICAL FILE

### M-CLASS CABINET



C/ Lanzahíta, nº 6-Pol.Ind. Albarreja  
28946 Fuenlabrada (Madrid)  
CIF: A-78025293  
Phone. 91.208.94.00  
Fax: 91.208.94.01

## GENERAL FEATURES:

- Cabinet with overbed table attached to the base. Structure in aluminium and HPL board.

- Sliding drawer in compact laminate and front in HPL



- Lower cupboard with 95° opening door. Front board made of HPL.

- Steel base including four swivelling Ø50 mm wheels, two of them with brakes.

- Chromed steel towel holder in the lateral side of the block.

- **OVERBED TABLE**

Aluminium structure and HPL board, tilting table in multiple positions. Manual height range adjustment from 79 to 108 cm. Folding overbed table along the side of the cabinet to reduce space when it is not used.



## TECHNICAL DATA:

OVERALL DIMENSIONS (W x D x H)	480 x 460 x 840 mm
TRAY	583 x 364 mm
DRAWER MEASURES	378 x 150 mm
CUPBOARD MEASURES	378 X 600 mm
OVERBED TABLE HEIGHT ADJUSTMENT	790 - 1080 mm
APPROXIMATE WEIGHT	30 KG.