

Width
70cm



46cm to 95cm



CLINNOVA Therapy 2 Section – Basic



Designed for head and neck mobilisation, as well as manual traction techniques, these strong and rigid tables feature a short head section, with the option of adjustable height armrests on the "Plus" model, providing additional patient comfort during prone treatments. The head section has a breathing hole and plug fitted as standard.

Models

See Page 2 for all available models



Model NV2558-CLS
Colour dark grey

Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Choice of hydraulic or electric height elevation
- Gas-assisted backrest adjustment from -25° to +45°
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 46cm to 95cm
- Linked retractable wheel design with adjustable foot
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

Optional Accessories

- 6200** Single Foot Switch
- 9002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 9005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 9001** Electric Perimeter Foot Switch (EPFS)



Negative head section to -25°



Offset gas strut for unobstructed patient view



Low patient transfer height

Models Available



NV2558-CLS

Hydraulic Elevation
Classic Base (CLS)



NV2568-CLS

Electric Elevation
Hand or Foot Switch
Classic Base (CLS)



NV2568-CLS-EPFS

Electric Elevation
Perimeter Foot Switch
Classic Base (CLS)



NV2558-PRM

Hydraulic Elevation
Premium Base (PRM)



NV2568-PRM

Electric Elevation
Hand or Foot Switch
Premium Base (PRM)



NV2568-PRM-EPFS

Electric Elevation
Perimeter Foot Switch
Premium Base (PRM)