



Medicare GP Gynaecology Couch



A multi-purpose couch providing electric positioning for gynaecology procedures, as well as a full length couch surface for most general prone and supine examinations and minor treatment procedures.

Models

- SM8543** Electric height, tilt and backrest with manual head section, retractable wheel system
- SM8543D** Deluxe Model complete with base cover, matching operators chair and head support cushion, retractable wheel system



Model SM8543 Colour black



Model 6045 optional gynaecology leg supports

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric adjustable angle backrest from horizontal to +80°
- Electrically adjustable seat angle with trendelenburg to -10°
- Power assisted adjustable head section from -25° to +50°
- Breathing hole and plug fitted
- Excellent height range from 50cm to 103cm
- Retractable wheels with adjustable foot for uneven surfaces
- Clearance beneath frame for mobile hoist
- Anti-bacterial powder coated frame
- Hand switch fitted as standard
- Supplied complete with lithotomy stirrups and debris tray
- Standard width upholstery 70cm
- Paper roll holder fitted as standard
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6053** Protective base frame cover
- 6057** Battery back-up facility
- 6045** Gynaecology leg supports
- 6007** Head support cushion
- 6209** Triple foot switch
- 6100** Matching operators chair

See page 48-49 for additional information on our range of optional accessories



See choice of matching operators chairs on page 50

Trendelenburg tilt to -10°. Model 6044 optional lithotomy stirrups