



## Sterling 3 Section Couch



Like the 2 Section Model, the Sterling 3 Section Couch is excellent value for money and the 3 section format provides a more versatile and comfortable position for seated procedures and examinations. The backrest support can be raised and fully adjusted from horizontal to +85° and the foot section is adjustable from -20° to +40° and can also be used as the head end for prone treatments with the option of adding a breathing hole and plug. The couch is therefore suitable for a wide range of examination and minor treatment procedures. A hand switch is fitted as standard for height control with a foot switch available as an option.

### Models

**SX3008** Electric with hand switch for height operation



Model SX3008 Sterling 3 section electric with optional breathing hole and plug. Colour black

### Standard Features

- Safe working load and lifting capacity of 200Kg
- Adjustable head/foot section from -20° to +40°
- Adjustable backrest from 0° to +85°
- Excellent height range from 46cm to 99cm
- 100mm diameter individual locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- Standard width upholstery - 63cm
- Choice of upholstery colours

### Optional Accessories

- 6005** Paper roll holder
- 6004** Breathing hole and plug
- 6003** Pair of side support rails
- 6003P** Pair of side support rail cushions
- 6200** Single foot switch operation
- 6007** Matching headrest/ knee support

See page 10 for additional information on our range of optional accessories.



Low patient transfer height