



ST4642



ST4642W



Therapy Deluxe Bobath Couch



The top of the range Deluxe Bobath couch has a safe working load and lifting capacity of 325Kg (715lbs). This model features electric height and backrest operation and has large central locking castors with steering facility for ease of movement. There are a wide range of specific accessories available to further enhance this model.

Models

- ST4642 Electric with electric backrest 105cm
- ST4642W Electric with electric backrest 125cm

Standard Features

- Safe working load and lifting capacity of 325Kg (715lbs)
- Electric height operation
- Electric backrest operation to +75°
- Large central locking castors with steering control
- Excellent height range from 46cm to 100cm
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Dual hand switch fitted as standard

Optional Accessories

- 6004 Breathing Hole & Plug
- 6057 Battery Back-Up
- 6239 Double Foot Switch
- 6001XL Perimeter Foot Switch
- 6074 Side Support Rails
- 6074P Side Support Rail Cushions
- 6041 Side Support Loops
- 6012 Upper and Lower Limb Support
- 6013 Patient Support Loops

See page 48-55 for additional information on our range of optional accessories and colour options



Model ST4642 Deluxe bobath 105cm, shown with optional side support rails and cushions. Colour lollipop



Model ST4642 with optional side support rails and hand support loops. Colour light grey 105cm wide



Model ST4642 with optional side support rails and patient limb support. Colour light grey 105cm wide



Model ST4642 with optional side support rails and hand support loops. Colour light grey 105cm wide



Model ST4642 with optional side support rails and patient limb support. Colour light grey 105cm wide