



## Therapy Mat Table



The Mat Therapy Tables are very similar to the Bobath models in terms of specification, but without the backrest facility. The single section top design provides a large and stable variable height support surface for a wide range of exercises and Bobath rehabilitation procedures.

### Models

<b>ST5551</b>	Hydraulic width 105cm
<b>ST5561</b>	Electric width 105cm
<b>ST5551W</b>	Hydraulic width 125cm
<b>ST5561W</b>	Electric width 125cm
<b>ST5561L</b>	Electric width 150cm



Model ST5561W Colour atomium



Model ST5561W Shown with optional side support rails and cushions. Colour atomium

### Standard Features

- Safe working load and lifting capacity of 250Kg (550lbs)
- Choice of widths available 105cm, 125cm or 150cm
- Choice of manual or electric height elevation on 105cm and 125cm models
- Excellent electric height range from 46cm to 100cm
- Retractable wheel design and adjustable foot for stability on uneven floors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- Bariatric models with a safe working load of 325Kg are available upon request

### Optional Accessories

- 6057** Battery Back-Up
- 6200** Single Foot Switch
- 6001XL** Perimeter Foot Switch
- 6074** Side Support Rails
- 6074P** Side Support Rail Cushions

See page 48-55 for additional information on our range of optional accessories and colour options



Model ST5551 Hydraulic mat table. Colour sky blue